# Instructions for Use Disposable Forceps (VS07XX.2X-XX)

**Product Description:** Vortex Surgical Disposable Forceps are a line of single use vitreoretinal forceps with various tip geometries. Each forceps has a handle with eight (8) individual levers allowing for 360 degree actuation and rotation. The forceps tip is opened and closed by depressing and releasing on the levers of the handle. The tail piece of the handle is removable depending on the surgeon's preference. The forceps are available in three gauges: 23ga (Orange), 25ga (Blue) and 27ga (Purple). All forceps are compatible with their corresponding gauge trocar cannula system.

**Indications for use:** Vortex Surgical Disposable Forceps are intended as a single-use instrument used in ophthalmic surgery for posterior segment procedures.

**Contraindications:** - None

**Known complications:** - Retinal detachment, tears, holes and contusions, bleeding, inflammation, and infection.

These complications are statistically rare and it is assumed that the user is adequately trained in the treatment of these known complications and methods of avoidance.

**Caution:** - Federal (USA) law restricts this device sale by, or on the order of, a physician.

**Instructions for use:** (Electronic Copy of IFU can be found at <https://www.vortexsurgical.com>)

- Determine package integrity.
- Remove product from protective packaging.
- Verify tip operation by actuating the handle and determine if tip is damaged. Do not attempt to repair or modify instrument.
- Affect sterile transfer of the product to the sterile field.
- Follow established surgical procedures.

**Device Storage:** Vortex Surgical, Inc recommends that the product is stored in a clean, dry and well-ventilated area at room temperature 15-37 °C (59-98.6 °F) away from direct sunlight.

**Sterilization:** This product is sterilized by Ethylene Oxide Gas (EtO) and provided with an EtO indicator. Verify indicator status color prior to use.

**Warnings:** Do not use if package integrity has been compromised.

Do not use if tip or handle is damaged.

Do not use if product exposed to conditions outside of indicated range.

Do not use if EtO indicator does not show exposure to sterilization gas.

**Reuse:** Do not reuse or reprocess instrument. Reuse could lead to cross contamination, infection, or device failure.

**Re-sterilization:** Do not re-sterilize instrument

**Device Disposal:** This single use surgical instrument should be considered clinical waste and should be disposed in accordance with clinical waste laws applicable in your country.

**Reporting:** Report serious incidents to Vortex Surgical and your regulatory body

Vortex Surgical, Inc. excludes all warranties, whether expressed, implied, or otherwise, to matters beyond the direct control of Vortex Surgical, Inc and the end result of this device's use. This would include, but not be limited to handling, shipment and storage of the device and patient diagnosis and treatment. The fitness and merchantability of this device are as specified. Implied factors are specifically excluded. Vortex Surgical, Inc. is not liable for loss, whether incidental or consequential, damage and/or expense, arising directly or indirectly from the use of this device. There is no additional liability or responsibility assumed, other than that specified. All additional liability or responsibility is specifically disclaimed.

 Vortex Surgical, Inc. 4 Research Park, Suite 124, Saint Charles, MO 63304, USA  
info@vortexsurgical.com 636-778-4350 www.vortexsurgical.com

Manufacturer


Sterilized using ethylene oxide    Keep dry    Do not reuse    Keep away from sunlight    Do not resterilize    Temperature limit 37°C

Do not use if damaged    Prescription use only    Authorized representative in the European Union    Consult Instructions for Use

Not made with natural rubber latex    Catalog number    Lot number    Quantity    Use By    Manufactured in USA/Date of Manufacture

Medical Device    2797    Vortex SURGICAL    © Vortex Surgical Inc., 2023

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